

Sales Conditions

Price Changes: Prices are subject to change without notice.

Currency:

Europe: Euro, EUR, excluding VAT.

Asia: Euro, EUR, excluding VAT.

USA: USD

Delivery Terms: EXW, Ex Works.

General Terms: ORGALIME S2000

Date of Delivery: The date given in our confirmations are only given in an indicative scope. Delays do not entitle client to cancel the order to refuse merchandise or to ask for damages and interest.

Payment Terms: 30 days net

VeruTech AB

Owner of the DiLab Brand
Galjevångsvägen 22
SE-224 65 Lund, Sweden
Tel: +46 46 211 2795
Fax: +46 46 37 37 11

Contact info:

E-mail: info@verutech.com
E-mail: info@dilab.com
www.dilab.com
www.verutech.com
twitter.com/verutech

AccuSampler[®]

Accessories,
Consumables &
Replacement parts

2012

VT
VERUTECH AB

DiLab[®]
www.dilab.com

Dear Customers,

Our main office is located in Sweden, and in the great majority of cases you should turn to us with your questions.



Galjevångsvägen 22
SE-224 65 Lund, Sweden
VAT-nr: SE556769648801
Phone: +46 (0) 46 211 2795
Fax: +46 (0) 46-373711
E-mail: info@verutech.com
<http://www.verutech.com>

To provide the best possible service for our customers we also have distributors for our AccuSampler® consumables in U.S and China. Information on how to contact these companies is provided below. Please note that all service and support goes through the VeruTech main office in Sweden.

Distributor (U.S)

Instech Solomon
Plymouth Meeting, PA USA
toll free 1-800-443-4227
tel +1-610-941-0132
fax +1-610-941-0134
email: orders@instechlabs.com






Distributor (China)

Shanghai Manpu Bio Technology Co.,LTD
Tel: +86-21-61361498
Fax: +86-21-61350931
[Http://www.manpubio.com.cn](http://www.manpubio.com.cn)
Add: No.250, Lane 399, Dushi Road, Minhang District, Shanghai, China





Table of Contents

Cages	1
Cage Accessories	2
Balance Arm	4
Swivels	5
Swivels & Accessories	7
Pumps & Accessories	9
Covance Harnesses & tethers.....	11
VAH™ for rats, single channel	13
VAH™ for rats, dual channel	15
VAH™ for mice, single channel	16
Buttons & tethers	17
Vascular Access Buttons™	18
Catheters for rats	19
Catheters for mice	20
Dual Tubings for AccuSampler® μ	21
Catheter accessories & connectors	22
Blunt needles	23
Clamps	24
Vials & Accessories	25
Fluids & Tube Systems	27
Training Materials	30
Wiring Diagrams.....	31
Sales Conditions	Back cover

Cages

Picture	Description	Article No.	Price
	<p>Animal cage - Large Used for rats. Includes: lid, bottom, locks, water bottle and holder.</p> <p><i>Material:</i> Polycarbonate <i>Amount:</i> 1 pcs</p>	100014	
	<p>Rodent Cage Black - Large For rodents like rats and guinea-pigs. Patent by Harlan Lab. <i>Production ended, available only when in stock.</i> <i>Material:</i> Polycarbonate <i>Amount:</i> 1 pcs</p>	100898	
	<p>Animal cage - Small Used for mice. Includes: lid, bottom, locks, water bottle and holder.</p> <p><i>Material:</i> Polycarbonate <i>Amount:</i> 1 pcs</p>	100188	
	<p>Rodent Cage Black - Small For small animals like mice. Includes: lid, bottom, locks, water bottle and holder</p> <p><i>Material:</i> Polycarbonate <i>Amount:</i> 1 pcs</p>	101087	
	<p>Complete metabolic cage Metabolic cage including stand and modified lid.</p> <p>For animals >300g</p> <p><i>Material:</i> Polycarbonate, Stainless steel <i>Amount:</i> 1 pcs</p>	100378	





Cage Accessories

Picture	Description	Article No. Price
	<p>Complete water bottle Standard water bottle with a Volume of 250 ml. Bottle and nipple.</p> <p><i>Material:</i> Polycarbonate, stainless steel, silicone gasket. <i>Amount:</i> 1 pcs</p>	100382
	<p>Holder for water bottle Fits Animal Cage Large or Small and Rodent Cage Black Large or Small. Includes 3 screws used to attach it to the cage.</p> <p><i>Material:</i> Hostaform, stainless steel <i>Amount:</i> 1 pcs</p>	100019
	<p>Metal cage lock Requires drilling in the cage (template/fixture included)</p> <p><i>Material:</i> Stainless Steel <i>Amount:</i> 2 pcs/pkg</p>	100922
	<p>Bottom for large animal cage The cage bottom separately. With metal cage locks premounted</p> <p><i>Material:</i> Polycarbonate & Stainless Steel <i>Amount:</i> 1 pcs</p>	100016
	<p>Top for large animal cage The cage top separately. With metal cage locks premounted</p> <p><i>Material:</i> Polycarbonate & Stainless Steel <i>Amount:</i> 1 pcs</p>	100686

Cage Accessories

Picture	Description	Article No.	Price
	<p>Complete lid for metabolic cage Simplifies the connection to the AccuSampler. Includes modified lid and lock.</p> <p><i>Material:</i> Polycarbonate, Stainless Steel and Hostaform <i>Amount:</i> 1 pcs</p>	DL-META-MOUNT	
	<p>Swivel mount for metabolic cages, slotted lid Slotted version for easier access to the animal.</p> <p><i>Material:</i> Hostaform <i>Amount:</i> 1 pcs</p>	DL-META-MOUNT/SL	
	<p>Urine Refrigeration Unit Used with the metabolic cage. For 30 ml bottles.</p> <p>For more information contact DiLab!</p>	100383	
	<p>Bottle for Urine Refrigerator Exchange bottle for Urine Refrigeration Unit.</p> <p><i>Size:</i> 30 ml bottles. <i>Material:</i> HDPE <i>Amount:</i> 1 pcs</p>	100606	

Balance Arm

Picture	Description	Article No.	Price
	<p>Balance arm Counter-balanced lever arm moves vertically and horizontally with the animal to prevent slack in the catheter (incl. Gimbal, supporting rod).</p> <p><i>Material:</i> Stainless Steel, aluminium <i>Amount:</i> 1 pcs</p>	DL-MCLA	
	<p>Mounting block for support rod Mounts the support rod on the animal cage.</p> <p><i>Material:</i> Stainless Steel & Aluminum <i>Amount:</i> 1 pcs</p>	DL-MCLA/MB	
	<p>Supporting rod collar This part is tightened around the supporting rod to keep the balance arm in position.</p> <p><i>Material:</i> Stainless Steel <i>Amount:</i> 1 pcs, minimum order 5 pcs</p>	DL-7060050028	
	<p>Supporting rod The supporting rod separately.</p> <p><i>Material:</i> Stainless Steel <i>Amount:</i> 1 pcs</p>	DL-7060050006	






Swivels

Picture	Description	Article No. Price
	<p>Dual channel swivel, 20 gauge Including 1 clamp (Article: DL-CLAMP). The second fluid channel gives the ability to infuse, withdraw and monitor simultaneously.</p> <p>To be used with the AccuSampler Standard instrument.</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/D/20
	<p>Dual channel swivel, 22 gauge Including 1 clamp (Article DL-CLAMP). The second fluid channel on these swivels gives you the ability to infuse, withdraw and monitor simultaneously.</p> <p>To be used with the AccuSampler μ instrument.</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/D/22
	<p>Single channel swivel, 25ga Including 1 plastic swivel to tether clamp. With standard inlet only. <i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/25PS
	<p>Single channel swivel, 22ga Including 1 plastic swivel to tether clamp. With standard or Luer inlet. <i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/22PS or DL-375/22PLS
	<p>Single channel swivel, 20ga Including 1 plastic swivel to tether clamp. With standard or Luer inlet. <i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/20PS or DL-375/20PLS

Picture	Description	Article No. Price
	<p>Single channel swivel, 18 gauge Including 1 clamp (article no. 100368).</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/18TW
	<p>Single channel swivel, 20 gauge Including 1 clamp (Article No 100368).</p> <p>Can be used with the AccuSampler Standard instrument.</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/20
	<p>Single channel swivel, 22 gauge Including 1 clamp (article no. 100368).</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/22
	<p>Single channel swivel, 25 gauge Including 1 clamp (article no. 100368).</p> <p><i>Material:</i> Stainless Steel, Teflon <i>Amount:</i> 1 pcs</p>	DL-375/25

Swivel Accessories

Picture	Description	Article No.	Price
	Swivel connector, 20 gauge Connects the yellow tubing with the swivel. Fits a 20 gauge swivel (Plastic or Stainless Steel). <i>Material:</i> PTFE <i>Amount:</i> 1 pcs	100267	
	Swivel connector, 22 gauge Connects the yellow tubing with the swivel. Fits a 22 gauge swivel (Plastic or Stainless Steel). <i>Material:</i> PTFE <i>Amount:</i> 1 pcs	100269	
	Swivel connector, luer Connects the yellow tubing with a swivel with luer inlet. <i>Material:</i> PTFE <i>Amount:</i> 1 pcs	100839	
	Rinsing Device For 20 gauge swivel. <i>Amount:</i> 1 pcs	100946	
	Rinsing Device For 22 gauge swivel. <i>Amount:</i> 1 pcs	100947	

Picture	Description	Article No. Price
	<p>Gimbal Holder for swivels. Fits both the 20 gauge (DL-375/D/20) and the 22 gauge swivel (DL-375/D/22).</p> <p><i>Material:</i> Hostaform, Stainless Steel <i>Amount:</i> 1 pcs</p>	DL-GIMBAL
	<p>Gimbal lock Fixation of gimbal to the balance arm (DL-MCLA).</p> <p><i>Material:</i> Hostaform <i>Amount:</i> 1 pcs</p>	DL-GIMBAL-LOCK
	<p>Gimbal wire retaining ring Fixation of steel wire on the gimbal.</p> <p><i>Material:</i> Hostaform <i>Amount:</i> 1 pcs</p>	DL-GIMBAL/WRR
	<p>Clamp, Universal To be used with the swivel. Connects the swivel to the tether. Mounts on the lower end of the swivel.</p> <p>Fits both the 20 gauge (DL-375/D20) and the 22 gauge swivel (DL-375/D/22).</p> <p><i>Material:</i> Stainless Steel, aluminium <i>Amount:</i> 5 pcs / pkg</p>	DL-CLAMP
	<p>Clamp, DiLab To be used with the swivel. Connecting the swivel to the tether. Mounts on the lower end of the swivel.</p> <p>Fits both the 20 gauge and the 22 gauge swivel.</p> <p><i>Material:</i> Hostaform <i>Amount:</i> 1 pcs</p>	100372

Pumps & Accessories


Picture	Description	Article No.	Price
---------	-------------	-------------	-------



Standard Sampler pump
Sample pump for AccuSampler®
Standard including connecting
cable.

100485

Calibrated and ready to use.



Pump head
Complete pump top for the
AccuSampler® Standard Samp-
ler Pump.

100486



Pressure Sensor
To be used with the Sampler pump
for the AccuSampler® Standard.
Uses advanced technique to
messure the pressure from the
animal.

100502







1 ml glass syringe
Approved to be used with the
AccuSampler μ (Service Pack 4)




101001

Amount: 1 pcs/pkg

Picture	Description	Article No. Price
	<p>Single syringe infusion pump 100374</p> <p><i>Features:</i></p> <ul style="list-style-type: none"> • High accuracy regarding flow and volume • Low flow rate • Even flow <p>Syringe size is specified at order: 1 and 5 ml are available.</p> <p>For AccuSampler Standard only</p>	
	<p>Dual syringe infusion pump 100373</p> <p>Injection device with 2 syringes for drug and saline delivery.</p> <p>After injection of the drug volume the saline syringe can immediately follow with a defined volume of saline. This ensures the fast injection of the drug into the animal and a safe flushing of the catheter.</p> <p>Syringe size is specified at order: 1 and 5 ml are available.</p> <p>For AccuSampler Standard only</p>	
	<p>1 ml Luer-Lok syringe 100535</p> <p><i>Amount:</i> 10 pcs/pkg</p>	
	<p>5 ml Luer-Lok syringe 100662</p> <p><i>Amount:</i> 10 pcs/pkg</p>	

Covance Harness & Tethers

Picture	Description	Article No. Price
	<p>Harness for rats, complete with spring for Single channel Including belly bands. Reusable and washable.</p> <p><i>Size:</i> Adjustable <i>Material:</i> Elastomer, Stainless steel <i>Amount:</i> 1 pcs</p>	DL-CIH95
	<p>Harness for rats, complete with spring for Dual channel Including belly bands. Reusable and washable.</p> <p><i>Size:</i> Adjustable <i>Material:</i> Elastomer, Stainless steel <i>Amount:</i> 1 pcs</p>	DL-CIH105
	<p>Rat Harness Including belly bands. Reusable and washable. Without tether. Single (95) or Dual channel (105) <i>Size:</i> Adjustable <i>Material:</i> Elastomer <i>Amount:</i> 1 pcs</p>	DL-CIH95AB or DL-CIH105AB
	<p>Bellybands for rat harness Seperate belly bands. Fits VAH95AB, CIH95 or CIH105 tethers.</p> <p><i>Size:</i> Adjustable <i>Material:</i> Elastomer <i>Amount:</i> 5 pcs/pkg</p>	DL-CIH95BV

Picture	Description	Article No. Price
	<p>Mouse Harness for Small Spring DL-CIH62 Including belly bands for single channel. Reusable and washable.</p> <p><i>Size:</i> Adjustable Spring ID:1,6mm <i>Material:</i> Elastomer & Stainless Steel <i>Amount:</i> 1 pcs</p>	
	<p>Mouse Harness for Large Spring DL-CIH62/PS95 Including belly bands for dual channels. Reusable and washable.</p> <p><i>Size:</i> Adjustable Spring ID:2,3mm <i>Material:</i> Elastomer & Stainless Steel <i>Amount:</i> 1 pcs</p>	
	<p>Bellybands for mouse harness DL-CIH62BV Seperate belly bands. 5 pcs/pkg.</p> <p><i>Size:</i> Adjustable <i>Material:</i> Elastomer <i>Amount:</i> 5 pcs/pkg</p>	


VAH™ for rats, Single channel

Picture	Description	Article No.	Price
	Vascular Access Harness for rats with catheter Including round tip catheter and SN22 needle <i>Size:</i> Adjustable <i>Material:</i> Silicone, Stainless steel, acetyl, elastomer <i>Amount:</i> 1 pcs	DL-VAH95AB-C	
	Vascular Access Harness for rats Including SN22 needle <i>Size:</i> Adjustable <i>Material:</i> Silicone, Stainless steel, acetyl <i>Amount:</i> 1 pcs	DL-VAH95AB	
	VAH rat tether Vascular Access harness tether <i>Size:</i> Adjustable <i>Material:</i> Polyurethane, stainless steel <i>Amount:</i> 1 pcs	DL-VAH95T	
	VAH tether kit, AccuSampler Std Including VAH tether and 20ga plastic swivel to fit the AccuSampler Standard. The kit is assembled and EtO sterilized. <i>Size:</i> Adjustable <i>Material:</i> Polyurethane, stainless steel <i>Amount:</i> 1 pcs/pkg	DL-KVAH95T/ DILAB	



Picture	Description	Article No.	Price
	Septum needle for VAH Extra rat VAH septum needles <i>Size: 22ga x 0,4in</i> <i>Amount: 12 pcs/pkg</i> or in bulk (50pcs/pkg)	DL-SN22 DL-SN22/ BULK50	
	Plug for VAH tether Sterile plug to seal off the tether connector and keep it clean while the animal is disconnected. <i>Material: Silicone, acetyl</i> <i>Amount: 10 pcs /pkg</i>	DL-SIP22P	
	Rat VAH Connector 3cm tubing, luer stub <i>Amount: 1 pcs</i>	DL-VAHLS22/30	

Instech's Vascular Access Harness™, an advancement on the Covance Infusion Harness™, permits quick, aseptic connection and disconnection of a catheterized rat and an infusion tether.



The system consists of two parts: a miniature external port housed in a harness (VAH95AB) and a mating spring tether with a connector (VAH95T) that pierces the septum to make the fluid connection.



VAH™ for rats, Dual channel

Picture	Description	Article No. Price
	<p>Vascular Access Harness for rats, two channel Including 2 catheters and SN22 needle</p> <p><i>Size:</i> Adjustable <i>Material:</i> Silicone, Stainless steel, elastomer, acetyl <i>Amount:</i> 1 pcs</p>	<p>DL-VAHD115-AB-C</p>
	<p>Vascular Access Harness for rats, two channel Including SN22 needle</p> <p><i>Size:</i> Adjustable <i>Material:</i> Silicone, Stainless steel, elastomer, acetyl <i>Amount:</i> 1 pcs</p>	<p>DL-VAH-D115AB</p>
	<p>VAH rat tether, two channel Dual channel Vascular Access harness tether assembly</p> <p><i>Size:</i> Adjustable <i>Material:</i> Polyurethane, stainless steel <i>Amount:</i> 1 pcs</p>	<p>DL-VAHD115T</p>
	<p>VAHD tether kit, AccuSampler Std Includes VAHD115T and a 20ga plastic swivel (375/20PS) to fit the AccuSampler Standard + a infusion port (SIP22/4) on the second channel. The kit is assembled and EtO sterilized. <i>Amount:</i> 1 kit/Pkg</p>	<p>KVAHD115T/DILAB</p>

VAH™ for mice, single channel

Picture	Description	Article No. Price
	<p>VAH harness for mice Injection set not included</p>	<p>DL-VAH62AB</p>
<p><i>Amount: 1 pcs</i></p>		
	<p>Mouse VAH injection set When accessing the harness port directly to infuse, withdraw or check patency, use this set rather than a septum needle for proper alignment and to avoid damage the small septum.</p>	<p>DL-VAHLS25/30</p>
<p>Mouse VAH tether assembly Vascular Access Harness Assembly for Mice. 7in</p>		
<p>DL-VAH62T</p>		
<p><i>Amount: 1 pcs</i></p>		
<p>Mouse VAH tether kit Vascular Access Harness Assembly for Mice + 25ga swivel, 24in CoEx and luer stub</p>		
<p>DL-KVAH62T</p>		
<p><i>Amount: 1 kit/pkg</i></p>		

Buttons & Tethers

Picture	Description	Article No.	Price
---------	-------------	-------------	-------



Dacron button

DL-DC105BS

Material: Dacron®
Amount: 10 pcs/pkg
ETO exposed



Dacron button

DL-DC105/
BULK

Material: Dacron®
Amount: 100 pcs/pkg



Spring Coil Tether

DL-PS105





Protects the catheters from the animal to the swivel.

Material: Stainless steel
Amount: 5 pcs / pkg
Length: 30 cm


Vascular Access Buttons™

Picture	Description	Article No. Price
	<p>Vascular Access Button for rats VAB with round-tip catheter and SN22 needle</p> <p><i>Material:</i> Polyester surgical felt, stainless steel, acetyl, silicone, polyurethane <i>Amount:</i> 1 pcs/pkg <i>ETO exposed</i></p>	VAB95BS-C
	<p>Vascular Access Button VAB with SN22 needle</p> <p><i>Material:</i> Polyester surgical felt, stainless steel, acetyl, silicone <i>Amount:</i> 1 pcs/pkg <i>ETO exposed</i></p>	VAB95BS
	<p>Vascular Access Button tether VAB tether assembly</p> <p><i>Material:</i> Polyurethane, stainless steel <i>Amount:</i> 1 pcs/pkg <i>ETO exposed</i></p>	VAB95T
	<p>Vascular Access Button kit VAB tether plus 20ga swivel to fit the AccuSampler Standard. The kit is assembled and EtO sterilized.</p> <p><i>Material:</i> Polyurethane, stainless steel <i>Amount:</i> 1 kit/pkg</p>	KVAB95T/ DILAB
	<p>Protective Cap for VAB Protective metal cap for Vascular Access Buttons.</p> <p><i>Material:</i> Aluminium <i>Amount:</i> 1 pcs/pkg</p>	VAB95CAP





Catheters for Rats

Picture	Description	Article No. Price
	<p>Rat Catheter for Carotid Artery & Jugular Vein With moveable beads.</p> <p>Production ended! This modell will be replaced with other combinations. <i>Contact us for more information</i> <i>ID: 0,38 mm OD: 0,76 mm</i> <i>Material: Tygon Length: 55 cm</i> <i>Amount: 10 pcs/pkg</i></p>	100936
	<p>Heparin coated catheter, Rat CBAS® for use with AccuSampler Standard. Rounded tip. Please note: Fits 22 gauge swivel without tubing connector or 20 gauge swivel together with tubing connector.</p> <p><i>Length: 60 cm</i> <i>OD: 0,9 mm ID:0,58mm</i> <i>Material: Polyurethane</i> <i>Amount: 1 pcs, minimum order 5 pcs</i> <i>ETO exposed</i></p>	DL-CBAS-C30
	<p>Silicon catheter for Rat Silicon is the old standard for long-term indwelling central venous catheters due to its softness and biocompatibility.</p> <p><i>Length: 60 cm</i> <i>OD: 0,94 mm ID:0,51mm</i> <i>Material: Silicone</i> <i>Amount: 1 pcs, minimum order 5 pcs</i> <i>ETO exposed</i></p>	DL-SIL-C30
	<p>Polyurethane catheter for Rat Polyurethan has supplanted silicon catheters due to its ease of insertation, durability and biocompatibility</p> <p><i>Length: 60 cm</i> <i>OD: 0,94 mm ID:0,58mm</i> <i>Material: PU</i> <i>Amount: 1 pcs, minimum order 5 pcs.</i> <i>ETO exposed</i></p>	DL-PU-C30




Catheters for Mice

Picture	Description	Article No. Price
	<p>Mouse Arterial Catheter With bead. Suitable with Dual Tubing no. 100766.</p> <p><i>Length:</i> 15,5 cm <i>Tip OD:</i> 0,36 mm <i>Material:</i> Polyurethane <i>Amount:</i> 1 pcs</p>	100828
	<p>Mouse Jugular Vein Catheter With bead. Suitable with Dual Tubing no. 100769.</p> <p><i>Length:</i> 9,5 cm <i>OD:</i> 0,76 mm <i>Material:</i> Polyurethane <i>Amount:</i> 1 pcs</p>	100583
	<p>Mouse Catheter for Jugular Vein and Carotid Artery With bead. Suitable with Dual Tubing no. 100769</p> <p><i>Length:</i> 10,5 cm <i>OD:</i> 0,6 mm <i>Material:</i> Silicon <i>Amount:</i> 1 pcs <i>ETO exposed</i></p>	100292
	<p>FunnelCath™ Mouse Catheter Tapers from 2 or 3Fr to 1,2Fr</p> <p><i>Total length:</i> 60 cm <i>Transition zone:</i> 4cm <i>1.2Fr zone:</i> 6cm <i>OD:</i> 0,6 mm <i>Material:</i> Polyurethane <i>Amount:</i> 1 pcs, minimum order 5 pcs. <i>ETO exposed</i></p>	DL-PUFC- C30-10

Dual Tubings for AccuSampler® μ

Picture	Description	Article No.	Price
	<p>Dual Tubing for AccuSampler® μ with 120 mm long common end Complete set from animal to swivel with a long common end to be implanted into the animal, as a catheter. For use on rats. <i>Patent by DiLab®.</i></p> <p><i>OD: 1 mm</i> <i>Material: Polyurethane</i> <i>Amount: 1 pcs</i> <i>ETO exposed</i></p>	100768	
	<p>Dual tubing for AccuSampler® μ with 120mm long common end Complete set from animal to swivel with a long common end to be implanted into the animal, as a catheter. For use on rats and mice. <i>Patent by DiLab®.</i> <i>OD: 0,84 mm</i> <i>Material: Polyurethane</i> <i>Amount: 1 pcs</i> <i>ETO exposed</i></p>	100287	
	<p>Dual Tubing for AccuSampler® μ with 25 gauge connector To be used together with catheter article no. 100828. <i>Patent by DiLab®.</i> <i>Material: Polyurethane</i> <i>Amount: 1 pcs</i> <i>ETO exposed</i></p>	100766	
	<p>Dual Tubing for AccuSampler® μ with 22 gauge connector To be used together with catheter article no. 100583 or 100292. <i>Patent by DiLab®.</i></p> <p><i>Material: Polyurethane</i> <i>Amount: 1 pcs</i> <i>ETO exposed</i></p>	100769	




Catheter, Accessories & connectors

Picture	Description	Article No.	Price
	Steel Connector, 22 gauge Non sterile <i>Length:</i> 8 mm <i>Material:</i> Stainless Steel <i>Amount:</i> 100 pcs/pkg	DL-SC22/8	
	Steel Connector, 25 gauge Non sterile <i>Length:</i> 10 mm <i>Material:</i> Stainless Steel <i>Amount:</i> 100 pcs/pkg	DL-SC25/10	
	Tubing Connector Connects a 22 ga swivel to suitable catheters. Do not need to be soaked in alcohol prior to installation. <i>Material:</i> Silicone <i>Length:</i> 10mm <i>ID:</i> 0,015in <i>Amount:</i> 100 pcs/pkg Sterile	DL-MC015/10	

Blunt Needles

Picture	Description	Article No.	Price
	Blunt Needle, 20 gauge Fits to catheter with article no. 100936. <i>Length: 1/2"</i> <i>Material:</i> Stainless Steel, polypropylene <i>Amount: 100 pcs/pkg</i>	DL-LS20	
	Blunt Needle, 22 gauge Fits to catheter with article no. 100583, 100292, Dual tubings with article no. 100287, 100766, 100769 and 100768. <i>Length: 1/2"</i> <i>Material:</i> Stainless Steel, polypropylene <i>Amount: 100 pcs/pkg</i>	DL-LS22	
	Blunt Needle, 25 gauge Fits to catheter with article no. 100828. <i>Length: 1/2"</i> <i>Material:</i> Stainless Steel, polypropylene <i>Amount: 100 pcs/pkg</i>	DL-LS25	
	Blunt Needle, 27 gauge Fits to catheter with article no. 100024 . <i>Length: 1"</i> <i>Material:</i> Stainless steel, polypropylene <i>Amount: 100 pcs/pkg</i>	DL-LS27	

Clamps

Picture	Description	Article No.	Price
	Vascular Clamp, Straight Bulldog clamps designed for vascular occlusion. <i>Material:</i> Radiopaque & Latex-Safe <i>Amount:</i> 2 pcs/pkg	100371	
	Vascular Clamp, Angled Bulldog clamps designed for vascular occlusion. <i>Material:</i> Radiopaque & Latex-Safe <i>Amount:</i> 2 pcs/pkg	100370	
	Catheter Clamp Catheter clip. <i>Length of jaw:</i> 4 mm <i>Material:</i> Stainless steel. <i>Amount:</i> 1 pcs	100387	

Vials & Accessories





Picture	Description	Article No.	Price
	Standard Vial Standard vial for AccuSampler Standard. Autoclaving max 20 min at 121 °C. Fits 16 and 32 position racks. <i>Max Volume:</i> 1,4 ml. <i>Material:</i> Polypropylene. <i>Amount:</i> 960 pcs/pkg	100384	
	Vial, Microtube for AccuSampler Standard Fits only 48 position racks. <i>Max Volume:</i> 300 µl. <i>Material:</i> PE <i>Amount:</i> 1000 pcs/pkg	100386	
	Microtainer™ LH Tubes Microtainer™ Tube, 200 - 400 µl, lithium heparin (no gel). Fits only 16 position racks. <i>Max Volume:</i> 400 µl <i>Amount:</i> 200 pcs/pkg	100385	
	Microtainer™ K2E Tubes Microtainer™ Tube, 250 - 500 µl, EDTA tube for whole blood, 0.8 mg EDTA (K2). Fits only 16 position racks. <i>Max Volume:</i> 500 µl <i>Amount:</i> 200 pcs/pkg	100642	
	Microtainer PST™ LH Tubes Microtainer™ tubes (with gel), 400-600µl. Fits only 16 position racks. <i>Max Volume:</i> 600 µl <i>Volume:</i> 200 pcs/pkg	100643	

Picture	Description	Article No. Price
	<p>Vial, 50μl Fits only AccuSampler μ.</p> <p><i>Max Volume:</i> 50 μl <i>Material:</i> Clear glass <i>Volume:</i> 100 pcs/pkg</p>	100291
	<p>Removable μ vial tray 48pos. Fits only AccuSampler$^{\circ}$ μ (SP4). Holds 48 vials for the μ-instrument</p> <p><i>Volume:</i> 1 pcs/pkg (vials not included)</p>	101090
	<p>Centrifuge adaptor This is a rotor adaptor for the blood sampling vials used in the Accusampler$^{\circ}$ μ. The adaptor fits most microcentrifuge tubes (Eppendorf) and contains a vial that encapsulates the micro vial.</p> <p><i>Amount:</i> 12 pcs/pkg</p>	100968
	<p>Gel loading tip Used to draw samples from the Vial, Microtube (article no 100291).</p> <p><i>Amount:</i> 200 pcs/pkg</p>	100641
	<p>Septum Lid for AccuSampler$^{\circ}$ Std, 16 positions</p>	100377
	<p>Septum Lid for AccuSampler$^{\circ}$ Std, 32 positions</p> <p>This option is designed to protect the samples in the rack from evaporation.</p> <p>Contact us for different rack types.</p>	100712

Fluids & Tube System

Picture	Description	Article No.	Price
	Sterile Water <i>Volume:</i> 100 ml/bag <i>Amount:</i> 40 pcs/pkg	100035	
	Saline Solution <i>Volume:</i> 100 ml/bag <i>Amount:</i> 40 pcs/pkg	100034	
	DiClean Cleaning solution for the tubing system of the AccuSampler® μ and Standard. Ready-to-use. Routine for cleaning the AccuSampler is described in the manual. <i>Volume:</i> 150 ml/bag <i>Amount:</i> 10 pcs/pkg	100909	
	Complete Tubing Set for AccuSampler® Standard Complete tubing set for AccuSampler® Standard, includes all necessary tubes and one green needle for the Saline bag (article no: 100268) and a 20g swivel connector (article no: 100267). <i>Volume:</i> 51-250 μ l <i>Material:</i> PFA Teflon	100006	

Picture	Description	Article No. Price
	<p>Complete Tubing Set for AccuSampler® μ Complete tubing set for AccuSampler® μ, includes all necessary tubes and two green needles for the Saline bag (article no: 100268) and two 22g swivel connectors (article no: 100269).</p> <p><i>Volume:</i> 5-50μ <i>Material:</i> PFA Teflon</p>	100183
	<p>Black Tube Black standard tubing for AccuSampler® Standard.</p> <p><i>Standard length:</i> 1800 mm <i>Material:</i> PFA Teflon <i>Amount:</i> 1 pcs</p>	100709
	<p>Brown Tube Brown standard tubing for AccuSampler® Standard</p> <p><i>Standard length:</i> 700 mm <i>Material:</i> PFA Teflon <i>Amount:</i> 1 pcs</p>	100698
	<p>Red Tube Red standard tubing for AccuSampler® Standard.</p> <p><i>Standard Length:</i> 100 mm <i>Material:</i> PFA Teflon <i>Amount:</i> 1 pcs</p>	100710

Picture	Description	Article No. Price
	<p>Yellow Tube Yellow tubing for AccuSampler® Standard. Available in standard (100708) or short version (100772). <i>Standard version:</i> 400 mm <i>Short version:</i> 225 mm <i>Material:</i> PFA Teflon <i>Amount:</i> 1 pcs</p>	<p>100708 or 100772</p>
	<p>Green Needle Cannula for connection to the bag with Saline Solution. Modified by VeruTech AB. <i>Dimension:</i> Ø 0,8 * 40 mm <i>Material:</i> Stainless steel <i>Amount:</i> 1 pcs</p>	<p>100268</p>
	<p>Septum for Waste Silicone rubber septum for the waste position. For AccuSampler® Standard <i>Material:</i> Rubber <i>Amount:</i> 10 pcs/pkg</p>	<p>100399</p>
	<p>Septum for Waste Silicone rubber septum for the waste position. For AccuSampler® µ <i>Material:</i> Rubber <i>Amount:</i> 10 pcs/pkg</p>	<p>101399</p>

Training Materials

Picture	Description	Article No.	Price
---------	-------------	-------------	-------



DVD - Permanent Catheterization of the Jugular Vein and Carotis Artery in the Mouse 100948



Manual for Mouse Catheterization 100949

Surgery instruction booklet, 32 pages.

DiLab[®] Wiring Diagram

AccuSampler[®] Std. - Rat



Swivel Connector 20 ga.
Article nr: 100267

Recommended setup



Dual channel swivel, 20ga
Article nr: DL-375/D/20



Single channel swivel, 20ga
Article nr: DL-375/PS20




Singe or Dual channel
Swivel, 20 gauge



Rat Catheter for Carotid
Artery & Jugular Vein
Article nr: 100936



Spring Coil Tether
Article nr: 100032



Vascular Access Harness
Tether
Article nr: DL-VAH95



Dacron Button, sterile
Article nr: DL-DC105



Harness for Rat
Article nr: DL-CH105



Vascular Access Harness
Article nr: DL-VAH

DiLab[®] Wiring Diagram AccuSampler[®] μ - Rat



Swivel Connector 22 ga.
Article nr: 100269



Swivel 22 ga.
Article nr: 100290



Dual Tubing,
120 mm, long common end
Article nr: 100287



Spring Coil Tether
Article nr: 100032




Dacron Button, sterile
Article nr: 100063



Harness for Rat
Article nr: 100654


DiLab[®] Wiring Diagram AccuSampler[®] μ - Mouse



Swivel Connector 22 ga.
Article nr: 100269



Swivel 22 ga.
Article nr: 100290



Dual Tubing
120 mm, long common end
Article nr: 100287



Dual Tubing
25 ga. connector
Article nr: 100766




Dual Tubing
22 ga. connector
Article nr: 100769



Mouse Arterial Catheter
Article nr: 100828




Mouse Jugular Vein
Catheter
Article nr: 100583



Mouse Carotid Artery
Catheter
Article nr: 100292



Light Weight Spring Coil
Tether
Article nr: 100968



Modified Mouse Harness
Article nr: 100869

